

For Research Use Only

# NDRG2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83902-4-PBS

## Basic Information

Catalog Number: 83902-4-PBS	GenBank Accession Number: BC010458	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57447	CloneNo.: 241039F6
Source: Rabbit	UNIPROT ID: Q9UN36	
Isotype: IgG	Full Name: NDRG family member 2	
Immunogen Catalog Number: AG2629	Calculated MW: 371 aa, 39 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00828-3, NDRG2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83902-6-PBS. Detection antibody: 83902-4-PBS. Standard: Ag2629. Range: 0.313-40 ng/mL