

For Research Use Only

CXCL8/IL-8 Recombinant antibody, PBS Only (Capture)

Catalog Number:83903-2-PBS



Basic Information

Catalog Number: 83903-2-PBS	GenBank Accession Number: BC013615	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3576	CloneNo.: 240889F6
Source: Rabbit	ENSEMBL Gene ID: ENSG00000169429	
Isotype: IgG	UNIPROT ID: P10145	
	Full Name: interleukin 8	
	Calculated MW: 99 aa, 11 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00830-2, CXCL8/IL-8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83903-2-PBS. Detection antibody: 83903-1-PBS. Standard: Eg0152. Range: 0.625-40 ng/mL.