

For Research Use Only

# CXCL8/IL-8 Recombinant antibody, PBS Only (Capture)

Catalog Number:83903-3-PBS



## Basic Information

Catalog Number:

83903-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC013615

GeneID (NCBI):

3576

ENSEMBL Gene ID:

ENSG00000169429

UNIPROT ID:

P10145

Full Name:

interleukin 8

Calculated MW:

99 aa, 11 kDa

Purification Method:

Protein A purification

CloneNo.:

240889H5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00830-1, CXCL8/IL-8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83903-3-PBS. Detection antibody: 83903-1-PBS. Standard: Eg0152. Range: 0.625-40 ng/mL.