

For Research Use Only

PLIN1 Recombinant antibody, PBS Only (Detector)



Catalog Number:83905-2-PBS

Basic Information

| | | |
|--------------------------------------|---|--|
| Catalog Number: 83905-2-PBS | GenBank Accession Number: NM_001145311 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 5346 | CloneNo.: 240994H2 |
| Source: Rabbit | UNIPROT ID: O60240 | |
| Isotype: IgG | Full Name: perilipin | |
| Immunogen Catalog Number: AG26426 | Calculated MW: 56 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00840-1, PLIN1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83905-1-PBS. Detection antibody: 83905-2-PBS. Standard: Ag26426. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00840-3, PLIN1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83905-3-PBS. Detection antibody: 83905-2-PBS. Standard: Ag26426. Range: 0.156-20 ng/mL