

For Research Use Only

# FBP2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83910-1-PBS

## Basic Information

Catalog Number: 83910-1-PBS	GenBank Accession Number: BC113632	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8789	CloneNo.: 240879B5
Source: Rabbit	UNIPROT ID: O00757	
Isotype: IgG	Full Name: fructose-1,6-bisphosphatase 2	
Immunogen Catalog Number: AG18245	Calculated MW: 339 aa, 37 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00855-3, FBP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83910-1-PBS. Detection antibody: 83910-4-PBS. Standard: Ag18245. Range: 0.313-40 ng/mL