

For Research Use Only

LOX Recombinant antibody, PBS Only (Capture)



Catalog Number: 83911-2-PBS

Basic Information

Catalog Number: 83911-2-PBS	GenBank Accession Number: BC074872	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4015	CloneNo.: 240945G4
Source: Rabbit	UNIPROT ID: P28300	
Isotype: IgG	Full Name: lysyl oxidase	
Immunogen Catalog Number: AG12479	Calculated MW: 417 aa, 47 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00857-1, LOX Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83911-2-PBS. Detection antibody: 83911-1-PBS. Standard: Ag12479. Range: 1.56-100 ng/mL