

For Research Use Only

NPT1 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number:83913-2-PBS

Basic Information

Catalog Number: 83913-2-PBS	GenBank Accession Number: BC101745	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6568	CloneNo.: 241010D3
Source: Rabbit	UNIPROT ID: Q14916	
Isotype: IgG	Full Name: solute carrier family 17 (sodium phosphate), member 1	
Immunogen Catalog Number: AG12819	Calculated MW: 467 aa, 51 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00846-2, NPT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83913-2-PBS. Detection antibody: 83913-4-PBS. Standard: Ag12819. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00846-3, NPT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83913-3-PBS. Detection antibody: 83913-2-PBS. Standard: Ag12819. Range: 0.625-80 ng/mL