

For Research Use Only

CHRNG Recombinant antibody, PBS Only (Capture)



Catalog Number:83915-3-PBS

Basic Information

Catalog Number: 83915-3-PBS	GenBank Accession Number: BC111802	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1146	CloneNo.: 240986D10
Source: Rabbit	UNIPROT ID: P07510	
Isotype: IgG	Full Name: cholinergic receptor, nicotinic, gamma	
Immunogen Catalog Number: AG16489	Calculated MW: 465 aa, 52 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00853-1, CHRNG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83915-3-PBS. Detection antibody: 83915-2-PBS. Standard: Ag16489. Range: 0.625-40 ng/mL



Cytometric bead array standard curve of MP00853-3, CHRNG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83915-3-PBS. Detection antibody: 83915-4-PBS. Standard: Ag16489. Range: 1.25-80 ng/mL