

For Research Use Only

NUCB2/nesfatin-1 Recombinant antibody, PBS Only (Detector)

Catalog Number:83917-4-PBS



Basic Information

Catalog Number: 83917-4-PBS	GenBank Accession Number: NM_005013	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4925	CloneNo.: 241037E5
Source: Rabbit	UNIPROT ID: P80303	
Isotype: IgG	Full Name: nucleobindin 2	
Immunogen Catalog Number: AG25152	Calculated MW: 50 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00864-2, NUCB2/nesfatin-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83917-1-PBS. Detection antibody: 83917-4-PBS. Standard: Ag25152. Range: 1.25-40 ng/mL