

For Research Use Only

# VISTA Recombinant antibody, PBS Only (Detector)



Catalog Number:83919-2-PBS

## Basic Information

Catalog Number: 83919-2-PBS	GenBank Accession Number: BC020568	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 64115	CloneNo.: 241003A12
Source: Rabbit	UNIPROT ID: Q9H7M9	
Isotype: IgG	Full Name: chromosome 10 open reading frame 54	
	Calculated MW: 311 aa, 34 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00865-1, VISTA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83919-1-PBS. Detection antibody: 83919-2-PBS. Standard: Eg0017. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00865-2, VISTA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83919-3-PBS. Detection antibody: 83919-2-PBS. Standard: Eg0017. Range: 0.156-20 ng/mL