

For Research Use Only

FAM185A Recombinant antibody, PBS Only (Detector)



Catalog Number:83921-2-PBS

Basic Information

Catalog Number: 83921-2-PBS	GenBank Accession Number: BC029175	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 222234	CloneNo.: 240804E2
Source: Rabbit	UNIPROT ID: Q8NOU4	
Isotype: IgG	Full Name: family with sequence similarity 185, member A	
Immunogen Catalog Number: AG15041	Calculated MW: 393 aa, 42 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00824-1, FAM185A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83921-1-PBS. Detection antibody: 83921-2-PBS. Standard: Ag15041. Range: 3.125-100 ng/mL