

For Research Use Only

CXCR7 Recombinant antibody, PBS Only



Catalog Number: 83927-1-PBS

Basic Information

Catalog Number: 83927-1-PBS	GenBank Accession Number: BC036661	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57007	CloneNo.: 240842C9
Source: Rabbit	UNIPROT ID: P25106	
Isotype: IgG	Full Name: chemokine (C-X-C motif) receptor 7	
Immunogen Catalog Number: AG14247	Calculated MW: 362 aa, 41 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83927-1-RR (CXCR7 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83927-1-PBS in a different storage buffer formulation.