

For Research Use Only

TIMP2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83938-1-PBS

Basic Information

Catalog Number: 83938-1-PBS	GenBank Accession Number: BC052605	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7077	CloneNo.: 241029F6
Source: Rabbit	UNIPROT ID: P16035	
Isotype: IgG	Full Name: TIMP metalloproteinase inhibitor 2	
	Calculated MW: 220 aa, 24 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00923-1, TIMP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83938-1-PBS. Detection antibody: 83938-2-PBS. Standard: Eg0430. Range: 1.25-80 ng/mL