

For Research Use Only

RAB11FIP5 Recombinant antibody, PBS Only (Capture)



Catalog Number:83939-3-PBS

Basic Information

Catalog Number: 83939-3-PBS	GenBank Accession Number: BC035013	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 26056	CloneNo.: 241038G10
Source: Rabbit	UNIPROT ID: Q9BXF6	
Isotype: IgG	Full Name: RAB11 family interacting protein 5 (class I)	
Immunogen Catalog Number: AG6128	Calculated MW: 70 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00924-1, RAB11FIP5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83939-3-PBS. Detection antibody: 83939-1-PBS. Standard: Ag6128. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP00924-2, RAB11FIP5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83939-3-PBS. Detection antibody: 83939-2-PBS. Standard: Ag6128. Range: 0.781-100 ng/mL