

For Research Use Only

CXXC5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83940-1-PBS

Basic Information

Catalog Number: 83940-1-PBS	GenBank Accession Number: BC017439	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 51523	CloneNo.: 241034B4
Source: Rabbit	UNIPROT ID: Q7LFL8	
Isotype: IgG	Full Name: CXXC finger 5	
Immunogen Catalog Number: AG9733	Calculated MW: 322aa, 33 kDa; 227aa, 24 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00925-1, CXXC5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83940-3-PBS. Detection antibody: 83940-1-PBS. Standard: Ag9733. Range: 0.313-40 ng/mL