

For Research Use Only

Mouse LAG-3/CD223 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83941-1-PBS



Basic Information

| | | |
|---------------------------------------|---|---|
| Catalog Number: 83941-1-PBS | GenBank Accession Number: NM_008479.2 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 16768 | CloneNo.: 241007B6 |
| Source: Rabbit | UNIPROT ID: Q61790 | |
| Isotype: IgG | Full Name: lymphocyte-activation gene 3 | |
| | Calculated MW: 57kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00927-1, MOUSE LAG-3/CD223 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83941-3-PBS. Detection antibody: 83941-1-PBS. Standard: Eg1220. Range: 0.313-40 ng/mL.