

For Research Use Only

Mouse CD30 ligand /Tnfsf8 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83942-1-PBS



Basic Information

Catalog Number: 83942-1-PBS	GenBank Accession Number: NM_009403.3	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 21949	CloneNo.: 241012C4
Source: Rabbit	UNIPROT ID: P32972	
Isotype: IgG	Full Name: tumor necrosis factor (ligand) superfamily, member 8	
	Calculated MW: 27 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00928-1, MOUSE CD30 ligand / Tnfsf8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83942-1-PBS. Detection antibody: 83942-2-PBS. Standard: Eg0667. Range: 0.313-40 ng/mL