

For Research Use Only

TMEM43 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83943-2-PBS

Basic Information

Catalog Number: 83943-2-PBS	GenBank Accession Number: BC008054	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 79188	CloneNo.: 241023B8
Source: Rabbit	UNIPROT ID: Q9BTV4	
Isotype: IgG	Full Name: transmembrane protein 43	
Immunogen Catalog Number: AG8522	Calculated MW: 400 aa, 45 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00929-1, TMEM43 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83943-1-PBS. Detection antibody: 83943-2-PBS. Standard: Ag8522. Range: 1.563-100 ng/mL.