

For Research Use Only

TOB1 Recombinant antibody, PBS Only (Detector)



Catalog Number:83944-3-PBS

Basic Information

Catalog Number: 83944-3-PBS	GenBank Accession Number: BC031406	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10140	CloneNo.: 241043C11
Source: Rabbit	UNIPROT ID: P50616	
Isotype: IgG	Full Name: transducer of ERBB2, 1	
Immunogen Catalog Number: AG6714	Calculated MW: 38 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00930-1, TOB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83944-4-PBS. Detection antibody: 83944-3-PBS. Standard: Ag6714. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00930-2, TOB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83944-1-PBS. Detection antibody: 83944-3-PBS. Standard: Ag6714. Range: 0.156-20 ng/mL