

For Research Use Only

EOMES Recombinant antibody, PBS Only (Capture)



Catalog Number:83945-4-PBS

Basic Information

Catalog Number: 83945-4-PBS	GenBank Accession Number: NM_001278182	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8320	CloneNo.: 241036G3
Source: Rabbit	UNIPROT ID: O95936	
Isotype: IgG	Full Name: eomesodermin homolog (Xenopus laevis)	
Immunogen Catalog Number: AG27634	Calculated MW: 73 aa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00931-1, EOMES Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83945-4-PBS. Detection antibody: 83945-3-PBS. Standard: Ag27634. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00931-2, EOMES Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83945-4-PBS. Detection antibody: 83945-1-PBS. Standard: Ag27634. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00931-3, EOMES Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83945-4-PBS. Detection antibody: 83945-2-PBS. Standard: Ag27634. Range: 0.156-20 ng/mL