

For Research Use Only

FANCD2 Recombinant antibody, PBS Only (Capture)



Catalog Number:83946-2-PBS

Basic Information

Catalog Number: 83946-2-PBS	GenBank Accession Number: NM_001018115	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2177	CloneNo.: 241048D8
Source: Rabbit	UNIPROT ID: Q9BXW9	
Isotype: IgG	Full Name: Fanconi anemia, complementation group D2	
Immunogen Catalog Number: AG29460	Calculated MW: 164 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00932-1, FANCD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83946-2-PBS. Detection antibody: 83946-3-PBS. Standard: Ag29460. Range: 0.781-100 ng/mL.