

For Research Use Only

LTBP4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83948-3-PBS



Basic Information

Catalog Number: 83948-3-PBS	GenBank Accession Number: NM_001042544	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8425	CloneNo.: 241056G6
Source: Rabbit	UNIPROT ID: Q8N2S1	
Isotype: IgG	Full Name: latent transforming growth factor beta binding protein 4	
Immunogen Catalog Number: AG34270	Calculated MW: 173 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00912-1, LTBP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83948-2-PBS. Detection antibody: 83948-3-PBS. Standard: Ag34270. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00912-2, LTBP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83948-1-PBS. Detection antibody: 83948-3-PBS. Standard: Ag34270. Range: 0.313-40 ng/mL



Bi-layer interferometry (BLI) kinetic assay of 83948-3-PBS against Human LTBP4 was performed. The affinity constant is 87.3 pM.