

For Research Use Only

RAB5A Recombinant antibody, PBS Only (Detector)



Catalog Number: 83949-2-PBS

Basic Information

Catalog Number: 83949-2-PBS	GenBank Accession Number: BC001267	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5868	CloneNo.: 241041H1
Source: Rabbit	UNIPROT ID: P20339	
Isotype: IgG	Full Name: RAB5A, member RAS oncogene family	
Immunogen Catalog Number: AG2549	Calculated MW: 215 aa, 24 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00914-1, RAB5A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83949-1-PBS. Detection antibody: 83949-2-PBS. Standard: Ag2549. Range: 0.313-100 ng/mL.



Cytometric bead array standard curve of MP00914-2, RAB5A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83949-3-PBS. Detection antibody: 83949-2-PBS. Standard: Ag2549. Range: 0.313-100 ng/mL.