

For Research Use Only

QSOX1/Quiescin Q6 Recombinant antibody, PBS Only (Capture)

Catalog Number:83952-2-PBS



Basic Information

Catalog Number: 83952-2-PBS	GenBank Accession Number: BC017692	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5768	CloneNo.: 240946D5
Source: Rabbit	UNIPROT ID: O00391	
Isotype: IgG	Full Name: quiescin Q6 sulfhydryl oxidase 1	
Immunogen Catalog Number: AG3396	Calculated MW: 747 aa, 83 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00903-1, QSOX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83952-2-PBS. Detection antibody: 83952-1-PBS. Standard: Ag3396. Range: 1.563-100 ng/mL.