

For Research Use Only

ATP6V1B1 Recombinant antibody, PBS Only (Capture)



Catalog Number:83953-4-PBS

Basic Information

Catalog Number: 83953-4-PBS	GenBank Accession Number: BC063411	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 525	CloneNo.: 240955F4
Source: Rabbit	UNIPROT ID: P15313	
Isotype: IgG	Full Name: ATPase, H ⁺ transporting, lysosomal 56/58kDa, V1 subunit B1	
Immunogen Catalog Number: AG6332	Calculated MW: 57 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00904-2, ATP6V1B1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83953-4-PBS. Detection antibody: 83953-1-PBS. Standard: Ag6332. Range: 0.781-100 ng/mL.