

For Research Use Only

C10orf33 Recombinant antibody, PBS Only (Capture)



Catalog Number:83954-3-PBS

Basic Information

Catalog Number: 83954-3-PBS	GenBank Accession Number: BC006131	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 84795	CloneNo.: 240984G10
Source: Rabbit	UNIPROT ID: Q8N2H3	
Isotype: IgG	Full Name: chromosome 10 open reading frame 33	
Immunogen Catalog Number: AG14701	Calculated MW: 581 aa, 63 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00905-2, C10orf33 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83954-3-PBS. Detection antibody: 83954-2-PBS. Standard: Ag14701. Range: 0.625-80 ng/mL.