

For Research Use Only

TLE3 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number:83955-2-PBS

Basic Information

Catalog Number: 83955-2-PBS	GenBank Accession Number: BC015729	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7090	CloneNo.: 240680D11
Source: Rabbit	UNIPROT ID: Q04726	
Isotype: IgG	Full Name: transducin-like enhancer of split 3 (E(sp1) homolog, Drosophila)	
Immunogen Catalog Number: AG1930	Calculated MW: 772 aa, 83 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00909-2, TLE3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83955-2-PBS. Detection antibody: 83955-1-PBS. Standard: Ag1930. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00909-3, TLE3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83955-3-PBS. Detection antibody: 83955-2-PBS. Standard: Ag1930. Range: 0.313-20 ng/mL