

For Research Use Only

XBP1 Recombinant antibody, PBS Only (Capture)



Catalog Number:83959-3-PBS

Basic Information

Catalog Number: 83959-3-PBS	GenBank Accession Number: BC000938	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7494	CloneNo.: 241004F6
Source: Rabbit	UNIPROT ID: P17861	
Isotype: IgG	Full Name: X-box binding protein 1	
Immunogen Catalog Number: AG21703	Calculated MW: 261 aa, 29 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00917-2, XBP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83959-3-PBS. Detection antibody: 83959-4-PBS. Standard: Ag21703. Range: 0.156-20 ng/mL