

For Research Use Only

SQRDL Recombinant antibody, PBS Only (Capture)



Catalog Number: 83960-1-PBS

Basic Information

Catalog Number: 83960-1-PBS	GenBank Accession Number: BC016836	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 58472	CloneNo.: 241226C12
Source: Rabbit	UNIPROT ID: Q9Y6N5	
Isotype: IgG	Full Name: sulfide quinone reductase-like (yeast)	
Immunogen Catalog Number: AG10762	Calculated MW: 450 aa, 50 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00918-1, SQRDL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83960-1-PBS. Detection antibody: 83960-2-PBS. Standard: Ag10762. Range: 0.78-100 ng/mL