

For Research Use Only

p57Kip2 Recombinant antibody, PBS Only (Capture)

Catalog Number:83966-3-PBS



Basic Information

Catalog Number: 83966-3-PBS	GenBank Accession Number: BC067842	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1028	CloneNo.: 241025D3
Source: Rabbit	UNIPROT ID: P49918	
Isotype: IgG	Full Name: cyclin-dependent kinase inhibitor 1C (p57, Kip2)	
Immunogen Catalog Number: AG19849	Calculated MW: 316 aa, 32 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00884-1, p57Kip2/CDKN1C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83966-3-PBS. Detection antibody: 83966-1-PBS. Standard: Ag19849. Range: 0.625-80 ng/mL.