

For Research Use Only

p57Kip2 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number:83966-4-PBS



Basic Information

Catalog Number:

83966-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19849

GenBank Accession Number:

BC067842

GeneID (NCBI):

1028

UNIPROT ID:

P49918

Full Name:

cyclin-dependent kinase inhibitor 1C
(p57, Kip2)

Calculated MW:

316 aa, 32 kDa

Purification Method:

Protein A purification

CloneNo.:

241025G11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00884-2, p57Kip2/CDKN1C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83966-2-PBS. Detection antibody: 83966-4-PBS. Standard: Ag19849. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00884-3, p57Kip2/CDKN1C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83966-4-PBS. Detection antibody: 83966-1-PBS. Standard: Ag19849. Range: 0.625-40 ng/mL