

For Research Use Only

SLC38A2 Recombinant antibody, PBS Only (Detector)



Catalog Number:83968-2-PBS

Basic Information

Catalog Number: 83968-2-PBS	GenBank Accession Number: BC040342	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 54407	CloneNo.: 241059D8
Source: Rabbit	UNIPROT ID: Q96QD8	
Isotype: IgG	Full Name: solute carrier family 38, member 2	
Immunogen Catalog Number: AG23082		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00889-1, SLC38A2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83968-1-PBS. Detection antibody: 83968-2-PBS. Standard: Ag23082. Range: 0.313-40 ng/mL.