

For Research Use Only

PRMT3 Recombinant antibody

Catalog Number:83969-3-RR



Basic Information

Catalog Number: 83969-3-RR	GenBank Accession Number: BC037544	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 10196	CloneNo.: 241092F2
Source: Rabbit	UNIPROT ID: O60678	Recommended Dilutions: IF/ICC 1:125-1:500
Isotype: IgG	Full Name: protein arginine methyltransferase 3	
Immunogen Catalog Number: AG13466	Calculated MW: 531 aa, 60 kDa	

Applications

Tested Applications: IF/ICC, FC (Intra), ELISA	Positive Controls: IF/ICC : HeLa cells,
Species Specificity: human	

Background Information

PRMT3(protein arginine methyltransferase 3) regulates protein functions by introducing asymmetric dimethylation marks at the arginine residues in proteins.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using PRMT3 antibody (83969-3-RR, Clone: 241092F2) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2).



1x10⁶ A549 cells were intracellularly stained with 0.25 ug PRMT3 Recombinant antibody (83969-3-RR, Clone:241092F2) and APC-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)(red), or 0.25 ug Isotype Control (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).