

For Research Use Only

Mouse CD40/TNFRSF5 Recombinant antibody, PBS Only (Detector)



Catalog Number:83972-1-PBS

Basic Information

| | | |
|--------------------------------|--|--|
| Catalog Number: 83972-1-PBS | GenBank Accession Number: NM_011611.2 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 21939 | CloneNo.: 240874A9 |
| Source: Rabbit | UNIPROT ID: P27512-1 | |
| Isotype: IgG | Full Name: CD40 antigen | |
| | Calculated MW: 32kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00896-1, MOUSE CD40/TNFRSF5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83972-2-PBS. Detection antibody: 83972-1-PBS. Standard: Eg1118. Range: 0.313- 40 ng/mL