

For Research Use Only

# Mouse CD40/TNFRSF5 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83972-4-PBS

## Basic Information

Catalog Number:

83972-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_011611.2

GeneID (NCBI):

21939

UNIPROT ID:

P27512-1

Full Name:

CD40 antigen

Calculated MW:

32kDa

Purification Method:

Protein A purification

CloneNo.:

240874G8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00896-2, MOUSE CD40/TNFRSF5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83972-3-PBS. Detection antibody: 83972-4-PBS. Standard: Eg1118. Range: 0.313- 40 ng/mL