

For Research Use Only

RNGTT Recombinant antibody, PBS Only (Detector)



Catalog Number:83974-3-PBS

Basic Information

Catalog Number: 83974-3-PBS	GenBank Accession Number: BC019954	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8732	CloneNo.: 240927D11
Source: Rabbit	UNIPROT ID: O60942	
Isotype: IgG	Full Name: RNA guanylyltransferase and 5'-phosphatase	
Immunogen Catalog Number: AG3105	Calculated MW: 597 aa, 69 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00899-1, RNGTT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83974-1-PBS. Detection antibody: 83974-3-PBS. Standard: Ag3105. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00899-2, RNGTT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83974-2-PBS. Detection antibody: 83974-3-PBS. Standard: Ag3105. Range: 0.781-100 ng/mL.