

For Research Use Only

P2RY13 Recombinant antibody, PBS Only (Capture)



Catalog Number:83975-1-PBS

Basic Information

| | | |
|---|---|---|
| Catalog Number: 83975-1-PBS | GenBank Accession Number: BC041116 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 53829 | CloneNo.: 240951D9 |
| Source: Rabbit | UNIPROT ID: Q9BPV8 | |
| Isotype: IgG | Full Name: purinergic receptor P2Y, G-protein coupled, 13 | |
| Immunogen Catalog Number: AG14108 | Calculated MW: 354 aa, 41 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00900-1, P2RY13 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83975-1-PBS. Detection antibody: 83975-2-PBS. Standard: Ag14108. Range: 0.781-100 ng/mL.