

For Research Use Only

# TIMM50 Recombinant antibody, PBS Only (Capture)



Catalog Number:83978-1-PBS

## Basic Information

Catalog Number:

83978-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17574

GenBank Accession Number:

BC121147

GeneID (NCBI):

92609

UNIPROT ID:

Q3ZCQ8

Full Name:

translocase of inner mitochondrial membrane 50 homolog (S. cerevisiae)

Calculated MW:

456 aa, 50 kDa

Purification Method:

Protein A purification

CloneNo.:

241097C12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00873-1, TIMM50 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83978-1-PBS. Detection antibody: 83978-1-PBS. Standard: Ag17574. Range: 0.313-40 ng/mL