

For Research Use Only

Mouse Serpin E1/PAI-1 Recombinant antibody, PBS Only (Capture)

Catalog Number:83980-1-PBS



Basic Information

Catalog Number:

83980-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

18787

UNIPROT ID:

P22777

Full Name:

serine (or cysteine) peptidase inhibitor, clade E, member 1

Calculated MW:

45kDa

Purification Method:

Protein A purification

CloneNo.:

241130E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00875-2, MOUSE Serpin E1/PAI-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83980-1-PBS. Detection antibody: 83980-4-PBS. Standard: Eg1100. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00875-3, MOUSE Serpin E1/PAI-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83980-1-PBS. Detection antibody: 83980-2-PBS. Standard: Eg1100. Range: 0.313-40 ng/mL