

For Research Use Only

Securin Recombinant antibody, PBS Only (Capture)

Catalog Number: 83981-1-PBS



Basic Information

Catalog Number:

83981-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12670

GenBank Accession Number:

BC101834

GeneID (NCBI):

9232

UNIPROT ID:

O95997

Full Name:

pituitary tumor-transforming 1

Calculated MW:

202 aa, 22 kDa

Purification Method:

Protein A purification

CloneNo.:

241090C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00876-2, Securin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83981-1-PBS. Detection antibody: 83981-4-PBS. Standard: Ag12670. Range: 1.56-100 ng/mL