

For Research Use Only

AGXT2L1 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83984-4-PBS

Basic Information

Catalog Number: 83984-4-PBS	GenBank Accession Number: BC022526	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 64850	CloneNo.: 241119G4
Source: Rabbit	UNIPROT ID: Q8TBG4	
Isotype: IgG	Full Name: alanine-glyoxylate aminotransferase 2-like 1	
Immunogen Catalog Number: AG10584	Calculated MW: 499 aa, 56 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00885-2, AGXT2L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83984-2-PBS. Detection antibody: 83984-4-PBS. Standard: Ag10584. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP00885-3, AGXT2L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83984-4-PBS. Detection antibody: 83984-3-PBS. Standard: Ag10584. Range: 0.625-80 ng/mL.