

For Research Use Only

NFE2L1 Recombinant antibody, PBS Only (Detector)



Catalog Number:83987-3-PBS

Basic Information

Catalog Number: 83987-3-PBS	GenBank Accession Number: BC010623	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4779	CloneNo.: 241089H5
Source: Rabbit	UNIPROT ID: Q14494	
Isotype: IgG	Full Name: nuclear factor (erythroid-derived 2)-like 1	
Immunogen Catalog Number: AG0525	Calculated MW: 772 aa, 85 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00878-1, NFE2L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83987-1-PBS. Detection antibody: 83987-3-PBS. Standard: Ag0525. Range: 0.625-80 ng/mL