

For Research Use Only

CRMP2 Recombinant antibody, PBS Only (Capture)

Catalog Number:83990-1-PBS



Basic Information

Catalog Number: 83990-1-PBS	GenBank Accession Number: BC067109	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1808	CloneNo.: 241088A5
Source: Rabbit	UNIPROT ID: Q16555	
Isotype: IgG	Full Name: dihydropyrimidinase-like 2	
Immunogen Catalog Number: AG5987	Calculated MW: 62 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00886-1, DPYSL2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83990-1-PBS. Detection antibody: 83990-2-PBS. Standard: Ag5987. Range: 6.25-100 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83990-1-PBS against Human CRMP2 was performed. The affinity constant is below 1 pM.