

For Research Use Only

E-cadherin Recombinant antibody, PBS Only (Detector)



Catalog Number:83991-2-PBS

Basic Information

Catalog Number: 83991-2-PBS	GenBank Accession Number: NM_004360.5	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 999	CloneNo.: 241098D10
Source: Rabbit	UNIPROT ID: P12830-1	
Isotype: IgG	Full Name: cadherin 1, type 1, E-cadherin (epithelial)	
	Calculated MW: 97 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00187-5, E-cadherin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83991-1-PBS. Detection antibody: 83991-2-PBS. Standard: Eg1112. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00187-6, E-cadherin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83991-4-PBS. Detection antibody: 83991-2-PBS. Standard: Eg1112. Range: 0.313-40 ng/mL.