

For Research Use Only

Histone H1.0 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number:83993-2-PBS

Basic Information

Catalog Number: 83993-2-PBS	GenBank Accession Number: BC000145	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3005	CloneNo.: 241147D11
Source: Rabbit	UNIPROT ID: P07305	
Isotype: IgG	Full Name: H1 histone family, member 0	
Immunogen Catalog Number: AG9982	Calculated MW: 21 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00895-2, Histone H1.0 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83993-3-PBS. Detection antibody: 83993-2-PBS. Standard: Ag9982. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP00895-3, Histone H1.0 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83993-2-PBS. Detection antibody: 83993-1-PBS. Standard: Ag9982. Range: 0.781-100 ng/mL.