

For Research Use Only

Mouse CTLA-4/CD152 Recombinant antibody, PBS Only (Detector)

Catalog Number:83995-4-PBS



Basic Information

Catalog Number: 83995-4-PBS	GenBank Accession Number: NM_009843.4	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 12477	CloneNo.: 241100G8
Source: Rabbit	UNIPROT ID: P09793	
Isotype: IgG	Full Name: cytotoxic T-lymphocyte-associated protein 4	
	Calculated MW: 25 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00907-1, MOUSE CTLA-4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83995-3-PBS. Detection antibody: 83995-4-PBS. Standard: Eg0632. Range: 0.156-20 ng/mL.