

For Research Use Only

NLRP7 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number:83996-2-PBS

Basic Information

Catalog Number: 83996-2-PBS	GenBank Accession Number: BC109125	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 199713	CloneNo.: 241120E4
Source: Rabbit	UNIPROT ID: Q8WX94	
Isotype: IgG	Full Name: NLR family, pyrin domain containing 7	
Immunogen Catalog Number: AG27465	Calculated MW: 112 kDa, 115 kDa, 118 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00908-1, NLRP7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83996-2-PBS. Detection antibody: 83996-3-PBS. Standard: Ag27465. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00908-3, NLRP7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83996-1-PBS. Detection antibody: 83996-2-PBS. Standard: Ag27465. Range: 0.313-40 ng/mL