

For Research Use Only

# iNOS Recombinant antibody, PBS Only (Detector)

Catalog Number:83997-2-PBS



## Basic Information

<b>Catalog Number:</b> 83997-2-PBS	<b>GenBank Accession Number:</b> BC130283	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 4843	<b>CloneNo.:</b> 230496H7
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P35228	
<b>Isotype:</b> IgG	<b>Full Name:</b> nitric oxide synthase 2, inducible	
<b>Immunogen Catalog Number:</b> AG17696	<b>Calculated MW:</b> 1153 aa, 131 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP00898-1, Human NOS2 Recombinant Matched Antibody Pair - PBS only. 83997-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag17696. 83997-2-PBS was HRP conjugated as the detection antibody. Range: 0.625-40 ng/mL.



Serum of fifteen individual healthy human donors was measured. The NOS2 concentration of detected samples was determined to be 9.62 ng/mL with a range of 0.10-94.17 ng/mL.