

For Research Use Only

CALCOCO1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 84009-3-PBS

Basic Information

Catalog Number:

84009-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13876

GenBank Accession Number:

BC003177

GeneID (NCBI):

57658

UNIPROT ID:

Q9P1Z2

Full Name:

calcium binding and coiled-coil
domain 1

Calculated MW:

691 aa, 77 kDa

Purification Method:

Protein A purification

CloneNo.:

241127G4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00945-1, CALCOCO1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84009-2-PBS. Detection antibody: 84009-3-PBS. Standard: Ag13876. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00945-2, CALCOCO1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84009-1-PBS. Detection antibody: 84009-3-PBS. Standard: Ag13876. Range: 0.313-40 ng/mL.