

For Research Use Only

# CALCOCO1 Recombinant antibody, PBS Only

Catalog Number:84009-4-PBS



## Basic Information

<b>Catalog Number:</b> 84009-4-PBS	<b>GenBank Accession Number:</b> BC003177	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 57658	<b>CloneNo.:</b> 241127G7
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9P1Z2	
<b>Isotype:</b> IgG	<b>Full Name:</b> calcium binding and coiled-coil domain 1	
<b>Immunogen Catalog Number:</b> AG13876	<b>Calculated MW:</b> 691 aa, 77 kDa	
	<b>Observed MW:</b> 100 kDa	

## Applications

**Tested Applications:**  
WB, IHC, Indirect ELISA

**Species Specificity:**  
human, mouse, rat

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)




W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84009-4-RR (CALCOCO1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84009-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84009-4-RR against Human CALCOCO1 were performed. The affinity constant is below 1 pM.